



Motion Computing Software Release Note

Name	Intel Mobile Clinical Assistant Win7
File Name	MCAPD_Setup.zip
Release Date	May 19, 2010
This Version	3.0.3.34538
Replaces Version	All previous versions
Intended Users	Motion C5v/F5v Users
Intended Product	Motion C5v/F5v
Supported Operating System	<ul style="list-style-type: none">• Windows XP Tablet Edition 2005*, Version 2002, Service Pack 2 or Service Pack 3 - Installed.• Windows 7 (32-bit)• * Tablet Edition 2005 is included with the free Service Pack 2 download. If you installed Windows XP Service Pack 2 on a Tablet PC, you now have Tablet Edition 2005.

Fixes in this Release

- Original Factory Release

What's New In This Release

- None

Known Issues and Limitations

- None

Installation and Configuration Instructions

To install this version update, perform the following steps:

1. Download **MCAPD_Setup.zip**
2. Download the Intel MCA 3.0.3.34538 to the directory location of your choice.
3. Unzip the file to a location on your local drive.
4. Double-click on the set-up.exe installer and follow the prompts.

Licensing

Refer to your applicable end-user licensing agreement for usage terms.